

APPLICATION NOTE

Analysis of 58 Volatile Organic Compounds Using Two CMS5000 Monitoring Systems

OVERVIEW

This application note describes the benefits of using the fully automated on-site CMS5000 Monitoring System to identify and alert users to the presence of volatile organic compounds (VOCs) in water systems. Fifty-eight VOCs are analyzed using two CMS5000 systems in tandem.

INTRODUCTION

On-site monitoring of VOCs is essential for the timely evaluation of potential hazards throughout water system infrastructures. VOCs are often inadvertently introduced into drinking water as the result of spills, improper disposal or soil migration¹. Many VOCs are known to cause severe health effects at low ppb concentrations².

CMS5000 employs SituProbe™ purge and trap sampling technology, gas chromatographic separation and Micro Argon Ionization Detection (MAID) to reliably monitor concentrations of multiple VOCs simultaneously in water, while in the presence of common interferences. The unit is designed to operate autonomously, analyzing numerous samples per day over several months without user supervision. The software generates an alarm if the concentrations exceed a user-defined threshold. Data collected from the calibrated CMS5000 allows users to make informed decisions based on their drinking water quality.

USEPA Methods 8260B and 5030 describe the analysis of approximately 90 VOCs in water by purge and trap. In order to analyze this number of compounds, the sample usually has to be analyzed using a mass spectrometer (MS). Most gas chromatography (GC) systems analyze 5-20 analytes. The separation of coeluting compounds when analyzing by GC is an analytical challenge.

In this application, 58 VOCs were analyzed on two separate CMS5000 systems with each employing complementary GC columns to improve separation of

co-eluting compounds. The systems were nearly identical, except one system used a non-polar 100% polydimethylsiloxane (PDMS) column and the other used a slightly polar 4% cyanopropyl-phenyl column.

EXPERIMENTAL

Calibration standards were prepared at 1, 3, 5 and 10 ppb by spiking 2 L of VOC-free water with the appropriate amount of a 200 µg/mL, 57 component VOC specialty mix in methanol and a 2000 µg/mL MTBE standard in methanol (SPEXCertiPrep). The standards were created by injecting 10.0 µL of the VOC specialty mix and 1.0 µL of the MTBE standard into 2 L of water and setting the concentrator fill time to 60, 180, 300 or 600 seconds.

Figure 1 and Figure 2 show chromatograms of the 1 ppb calibration standard. Using argon carrier gas, the analytes were purged from the calibration standard into the sampling tube headspace. The headspace sample was collected onto the Tri-Bed Concentrator. The analytes were thermally desorbed from the concentrator and then separated on a capillary column using a 51 minute temperature programmed method. A four point calibration curve with quadratic fit was generated from the data.

Calibration accuracy was measured by analyzing a 2 ppb calibration verification standard and calculating the percent recovery for each of the 58 compounds. See Table 1 for a list of compounds, the percent recovery, and the retention time of the compounds.

NOTE: MTBE, 2,2-dichloropropane, n-propylbenzene, 1,1,1,2-tetrachloroethane and 1,1,2,2-tetrachloroethane coeluted with a different compound on each column and could not be quantified by the software. The concentration of these compounds can be manually calculated by subtracting the concentration of the coeluting compound, determined on one column, from the concentration of the coeluting compounds on the other column. See Table 2.

CONCLUSION

The fully-automated CMS5000 Monitoring System is ideally suited to continually monitor VOCs in drinking water at low concentrations. The calibrated CMS5000 successfully separated, identified and quantified 58 VOCs at 1 to 10 ppb. Running an analysis with two complementary GC columns, with different stationary phases, reduced the number of coelutions and provided quantifiable data. These analyses demonstrate the utility of the CMS5000 for alerting users to unacceptable VOC concentrations in water, allowing water utilities to make

fast, informed decisions to protect and ensure public health and well-being.

REFERENCES

- 1 <http://www.health.state.mn.us/divs/eh/hazardous/topics/vocs.html>
- 2 Zogorski, J. S. (2006). Volatile Organic Compounds in the Nation's Ground Water and Drinking-Water Supply Wells. Reston, VA: U.S. Department of the Interior U.S. Geological Survey.

Figure 1 Chromatogram of 58 compounds at 1.0 ppb run on a DB-1 column




Figure 2 Chromatogram of 58 compounds at 1.0 ppb run on a DB-624 column




Table 1 Recovery and retention times of 58 compound standard

| No. | Compound | RT DB-1 | % Rec DB-1 | RT DB-624 | % Rec DB-624 |
|-----|---------------------------|---------|------------|-----------|--------------|
| 1 | vinyl chloride | 2:27 | 93 | 1:48 | NA |
| 2 | bromomethane | 2:55 | 168 | 1:55 | NA |
| 3 | chloroethane | 3:07 | 96 | 2:08 | 99 |
| 4 | trichlorofluoromethane | 3:58 | 125 | 2:20 | 97 |
| 5 | 1,1-dichloroethene | 4:46 | 92 | 2:47 | 91 |
| 6 | methylene chloride | 4:46 | 92 | 3:17 | 97 |
| 7 | trans-1,2-dichloroethene | 6:17 | 99 | 3:34 | 93 |
| 8 | 1,1-dichloroethane | 6:38 | 97 | 4:07 | 94 |
| 9 | MTBE | 6:38 | 97 | 3:34 | 93 |
| 10 | cis-1,2-dichloroethene | 8:12 | 100 | 5:01 | 82 |
| 11 | bromochloromethane | 8:33 | 103 | 5:29 | 115 |
| 12 | chloroform | 8:52 | 93 | 5:39 | 104 |
| 13 | 2,2-dichloropropane | 8:52 | 93 | 5:01 | 82 |
| 14 | 1,2-dichloroethane | 10:42 | 99 | 6:48 | 100 |
| 15 | 1,1,1-trichloroethane | 11:26 | 100 | 5:56 | 103 |
| 16 | 1,1-dichloropropene | 12:26 | 101 | 6:18 | 93 |
| 17 | benzene | 12:54 | 101 | 6:48 | 100 |
| 18 | carbon tetrachloride | 13:25 | 96 | 6:18 | 93 |
| 19 | dibromomethane | 15:59 | 99 | 10:00 | 99 |
| 20 | 1,2-dichloropropane | 15:59 | 99 | 9:32 | 103 |
| 21 | bromodichloromethane | 16:55 | 103 | 10:50 | 100 |
| 22 | trichloroethene | 16:55 | 103 | 8:42 | 104 |
| 23 | cis-1,3-dichloropropene | 19:44 | 103 | 13:07 | 103 |
| 24 | trans-1,3-dichloropropene | 21:15 | 98 | 16:14 | 102 |
| 25 | 1,1,2-trichloroethane | 21:40 | 102 | 16:14 | 102 |
| 26 | toluene | 22:26 | 101 | 14:38 | 101 |
| 27 | 1,3-dichloropropane | 22:26 | 101 | 17:49 | 103 |
| 28 | dibromochloromethane | 23:29 | 91 | 18:52 | 103 |
| 29 | 1,2-dibromoethane | 24:06 | 102 | 19:08 | 103 |
| 30 | tetrachloroethene | 25:21 | 100 | 17:17 | 98 |
| 31 | 1,1,1,2-tetrachloroethane | 27:02 | 103 | 22:30 | 101 |
| 32 | chlorobenzene | 27:02 | 103 | 21:43 | 103 |
| 33 | ethylbenzene | 28:02 | 100 | 22:30 | 101 |
| 34 | m-xylene | 28:30 | 95 | 23:09 | 101 |
| 35 | p-xylene | 28:30 | 95 | 23:09 | 101 |
| 36 | bromoform | 28:30 | 95 | 25:58 | 99 |
| 37 | styrene | 29:23 | 92 | 25:09 | 109 |
| 38 | o-xylene | 29:39 | 92 | 24:57 | 105 |
| 39 | 1,1,2,2-tetrachloroethane | 29:39 | 92 | 29:01 | 87 |
| 40 | 1,2,3-trichloropropane | 29:58 | 100 | 29:01 | 87 |

Table 1 Recovery and retention times of 58 compound standard (continued)

| No. | Compound | RT DB-1 | % Rec DB-1 | RT DB-624 | % Rec DB-624 |
|-----|------------------------|---------|------------|-----------|--------------|
| 41 | isopropylbenzene | 31:10 | 100 | 27:02 | 98 |
| 42 | bromobenzene | 31:20 | 103 | 28:04 | 105 |
| 43 | 2-chlorotoluene | 32:23 | 100 | 29:01 | 87 |
| 44 | 4-chlorotoluene | 32:33 | 98 | 29:41 | 103 |
| 45 | n-propylbenzene | 32:33 | 98 | 29:01 | 87 |
| 46 | 1,3,5-trimethylbenzene | 33:09 | 102 | 32:45 | 98 |
| 47 | 1,2,4-trimethylbenzene | 34:17 | 98 | 31:29 | 101 |
| 48 | tert-butylbenzene | 34:17 | 98 | 31:17 | 98 |
| 49 | 1,3-dichlorobenzene | 34:41 | 103 | 32:26 | 103 |
| 50 | 1,4-dichlorobenzene | 34:52 | 103 | 34:21 | 103 |
| 51 | sec-butylbenzene | 35:03 | 102 | 32:15 | 95 |
| 52 | 4-isopropyltoluene | 35:32 | 99 | 33:51 | 90 |
| 53 | 1,2-dichlorobenzene | 35:51 | 96 | 34:21 | 103 |
| 54 | n-butylbenzene | 36:49 | 102 | 34:44 | 96 |
| 55 | 1,2,4-trichlorobenzene | 41:23 | 103 | 41:00 | 102 |
| 56 | naphthalene | 41:44 | 102 | 41:50 | 97 |
| 57 | 1,2,3-trichlorobenzene | 42:42 | 105 | 42:51 | 100 |
| 58 | hexachlorobutadiene | 43:05 | 102 | 41:50 | 97 |

Table 2 Determination of the concentration of compounds which coelute on both columns

| Conc. of Cmpd. (ppb) | Equation |
|-------------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| MTBE (1) | (Conc. of MTBE + 1,1-dichloroethane) DB-1 - (Conc. of 1,1-dichloroethane) DB-624 = Conc of MTBE |
| MTBE (2) | (Conc. of MTBE + trans-1,2-dichloroethene) DB-624 - (Conc. of trans-1,2-dichloroethene) DB-1 = Conc. of MTBE |
| 2,2-dichloropropane (1) | (Conc. of 2,2-dichloropropane + chloroform) DB-1 - (Conc. chloroform) DB-624 = Conc of 2,2-dichloropropane |
| 2,2-dichloropropane (2) | (Conc. of 2,2-dichloropropane + cis-1,2-dichloroethene) DB-624 = Conc. of 2,2-dichloropropane |
| 1,1,1,2-tetrachloroethane (1) | (Conc. of 1,1,1,2-tetrachloroethane + chlorobenzene) DB-1 - (Conc. of chlorobenzene) DB-624 = Conc. of 1,1,1,2-tetrachloroethane |
| 1,1,1,2-tetrachloroethane (2) | (Conc. of 1,1,1,2-tetrachloroethane + ethylbenzene) DB-624 - (Conc. of ethylbenzene) DB-1 = Conc. of 1,1,1,2-tetrachloroethane |
| 1,1,2,2-tetrachloroethane | (Conc. of 1,1,2,2-tetrachloroethane + o-xylene) DB-1 - (Conc. of o-xylene) DB-624 = Conc. of 1,1,2,2-tetrachloroethane |
| n-propylbenzene | (Conc. of n-propylbenzene + 4-chlorotoluene) DB-1 - (Conc. of 4-chlorotoluene) DB-624 = Conc. of n-propylbenzene |



www.inficon.com reachus@inficon.com

Due to our continuing program of product improvements, specifications are subject to change without notice. All trademarks are the property of their respective owners.

diag53a1 ©2017 INFICON