

Micro GC Fusion® - Monitoring Hydrogen Concentrations in Ambient Air

Introduction

This datasheet highlights the 10 m Molsieve module used to test for hydrogen in ambient air. Several chromatograms are overlaid below - each with a different concentration of hydrogen.

Starting Parameters

These parameters can be used as a starting point for creating a method and can be adjusted to ensure all compounds are fully separated. Exact retention times will vary from GC to GC, but the compound order remains the same.

| Method Parameter | Module – 10 m Molsieve, Backflush Injector (GCM-W02) |
|---------------------------------|--|
| Inject time | 60 ms |
| Backflush time | 60 s |
| Injector temperature | 90°C |
| TCD temperature | 70°C |
| Column pressure and carrier gas | 40 psi, 99.999% nitrogen |
| Column temperature | 110°C (30 s) |
| Sample pump time | 15 s |
| Sample inlet temperature | 50°C |

Chromatogram

