

Research Crystals

MAXIMUM STABILITY IN DEMANDING APPLICATIONS

The standard AT-cut, 1 in. (25.4 mm) diameter crystal, pioneered by INFICON, is used primarily in research applications with the INFICON RQCM (Research Quartz Crystal Microbalance). These research crystals are designed for use with INFICON's rugged, easy to use crystal holders and are used for in-situ monitoring of deposited films or liquids.

Geometry: AT-cut quartz has superior mechanical and piezoelectric properties. The 1 in. (25.4 mm) diameter allows optimal distance between the active area of the crystal and the mounting o-ring. This improves the overall stability of the crystal by reducing the frequency changes due to mounting stress.

Finish: INFICON 1 in. (25.4 mm) research crystals are available with either an unpolished or polished finish. Polished crystals allow for a good correlation between theory and measurement during liquid immersion experiments.

APPLICATIONS

- Metal plating
- Gas-phase mass detection
- Electro-chemistry
- Micro-rheology
- Protein adsorption/desorption
- Biofilm formation

FEATURES AT A GLANCE

- Accurate
- Reliable
- Stable
- Wide variety of electrode materials available

Electrode Material: INFICON 1 in. (25.4 mm) research crystals are available in a variety of electrode materials including Gold, Platinum, and Titanium. Crystals with custom electrode materials are available on request. INFICON also offers Gold electrode crystals with an additional SiO₂ outer layer to create a hydrophilic surface needed for some biological applications.

Temperature: Even though AT-cut crystals are designed to minimize the change in frequency due to temperature, the effect of temperature can be significant when attempting to resolve small mass (frequency) changes over long periods of time. INFICON offers AT-cut crystals optimized for two operating temperatures, 90°C and 25°C. These crystals have very good frequency stability when operating close to their specified temperature (see Figure 1).

Contact INFICON today to request a quote for custom crystals to accommodate your special requirements.

SPECIFICATIONS

Frequencies:	5 MHz and 9 MHz
Frequency Range:	4.976 – 5.020 MHz 8.976 – 9.036 MHz
Resistance:	5 MHz, ~10 Ohms 9 MHz, ~7 Ohms
Crystal Diameter:	1 in. (25.4 mm)
Crystal Contour:	Plano-Plano
Crystal Orientation:	AT-Cut, for frequency stability over temperature
Compatibility:	INFICON/Maxtek liquid probes and crystal holders
Front Electrode Diameter:	0.5 in. (12.7 mm)
Packaging:	Individually packaged/sold
Crystal Thickness:	5 MHz, 0.013 in. (333 μm) 9 MHz, 0.007 in. (185 μm)
Surface Roughness:	50 Å (polished)

Figure 1



ORDERING INFORMATION

INFICON 1 in. (25.4 mm) Diameter Crystals - Electrode Configuration





REAR SIDE (CONTACT ELECTRODE)

1 in. (25.4 mm) Crystals, AT-Cut, 90°C						
Part No./Finish		F	Electrode Material (Front)			
Polished	Unpolished	Frequency	Adhesion Layer	Sensing Electrode		
149211-1	149211-2	5 MHz	Chromium	Gold		
149238-1	149238-2	5 MHz	Titanium	Gold		
149240-1	149240-2	5 MHz	Titanium	Platinum		
149241-1	149241-2	5 MHz	Chromium	Platinum		
149242-1	_	5 MHz	Titanium	Silver		
149243-1	_	5 MHz	Titanium	Aluminum		
149248-1	—	5 MHz	Titanium (No adhesion layer)			
149250-1	_	5 MHz	Chromium	Nickel		
149252-1	—	5 MHz	Chromium	Iron		
149255-1	_	9 MHz	Chromium	Gold		
149270-1	—	5 MHz	Titanium	SiO ₂		
149280-1	-	5 MHz	Chromium (No adhesion layer)			
750-1030-G1	—	5 MHz	Aluminum (No adhesion layer)			
750-1031-G1	_	5 MHz	Uncoated (Chromium/Gold contact electrode and wrap-around)			
750-1035-G1	-	5 MHz	Uncoated (Titanium/Gold contact electrode and wrap-around)			
SPC-1193-P1	-	5 MHz	Titanium/Gold	Iron		

1 in. (25.4 mm) Crystals, AT-Cut, 25°C							
Part No./Finish		Francis	Electrode Material (Front)				
Polished	Unpolished	Frequency	Adhesion Layer	Sensing Electrode			
149257-1	_	5 MHz	Titanium	Gold			
149272-1	—	9 MHz	Chromium	Gold			
149273-1	_	5 MHz	Chromium	Gold			
149274-1	_	5 MHz	Titanium	Platinum			
149275-1	_	5 MHz	Titanium	Silver			
149276-1	_	5 MHz	Titanium	Aluminum			
149277-1	149277-2	5 MHz	Titanium/Gold/Titanium	SiO ₂			
750-1029-G1	750-1033-G1	5 MHz	Titanium (No adhesion layer)				
750-1037-G1	_	9 MHz	Titanium	Platinum			
750-1039-G1	_	9 MHz	Titanium/Gold/Titanium	SiO ₂			
149296-1	_	9 MHz	Titanium	Gold			



www.inficon.com reachus@inficon.com Due to our continuing program of product improvements, specifications are subject to change without notice. All trademarks are the property of their respective owners. cibf15a1-a ©2013 INFICON