

# EZ IQ Hardware and Software Requirements

## Installation Requirements for EZ IQ

### Introduction

EZ IQ is validated on both the English and Western European Language versions of Windows.<sup>®</sup>

### Computer Minimum Requirements

Type	Microsoft <sup>®</sup> Operating System	CPU	RAM <sup>1</sup>	Disk Space	Network	Graphics
<b>Minimum Configuration</b>	XP Pro. SP3 with .Net 3.0 or Vista Business SP1 or Windows 7 Pro 32-bit	2 GHz	2 GB	10 GB	See instrument requirements	1024 x 768
<b>Recommended Configuration</b>	XP Pro. SP3 with .Net 3.0 or 3.5 SP1 or Vista Business SP1 or Windows 7 Pro 32-bit	Dual Core 2 GHz	2 GB	10 GB and >60% free HD space	100 Mbps or Higher	1280 x 1024 or higher

<sup>1</sup> RAM total available memory for EZ IQ including the operating system, but no other application.

### Important Notes

- ◆ Microsoft TCP/IP networking is required for all products. WAN's (wide area networks) are not supported.
- ◆ Disk Space requirements should be adjusted based on archival periodicity. INFICON recommends a minimum of one year of the expected disk usage on top of the recommendation of the operating system. For reference, each analysis will generate a data file of ~200Kb.

### Upgrade Support

Micro GC methods and data defined and generated under Class Driver cannot be migrated to the Enhanced Driver. Both driver types can co-exist but not co-execute when controlling the same instrument. See online help for details.

